

Analysis of SNR in AMDX2011P retinal images

Alex Matov
November 2022

Types of pixels in the image & computation approach

Very bright pixels (type of hot pixels) – excluded

Protein aggregates (AMDX-2011P) – foreground





Cellular fluorescence (non-specific labeling) – background

Camera noise (dark noise, shot noise) – excluded

- **Areas of foreground & background in each image are identified by visual inspection (Julie)**
- **I segment the images based on this input into two intensity ranges (background & foreground)**
- **In control images, regions selected by eye (Julie) have intensity ranges of 50-100 for foreground & 10-50 for background (see histograms)**
- **In TTR images, regions selected by eye (Julie) have intensity ranges of 80-150 for foreground & 10-30 for background (see histograms)**
- **The average of the foreground pixels over the average of the background pixels is the SNR, which is 3-fold higher in TTR**





frgr = pixels with intensity between 80 & 150
bkgr = pixels with intensity between 10 & 30

SNR analysis of TTR images



frgr = pixels with intensity between 80 & 150
bkgr = pixels with intensity between 10 & 30

SNR analysis of TTR images



48 TTR ON

SNR analysis of TTR images

counts

Intensity [A.U]

frgr = pixels with intensity between 80 & 150

bkgr = pixels with intensity between 10 & 30



Background 1/0 mask – pixels with intensity 10-30

Foreground 1/0 mask – pixels with intensity 80-150


Raw ex vivo AMDX-2011 8-bit image by Julie

SNR analysis of TTR images

SNR = 5.92 (the average pixel intensity in the foreground over the average pixel intensity in the background)




Raw ex vivo AMDX-2011 8-bit image
by Julie




SNR analysis of TTR images


frgr = pixels with intensity between 80 & 150
bkgr = pixels with intensity between 10 & 30



Foreground 1/0 mask – pixels with intensity 80-150




Background 1/0 mask – pixels with intensity 10-30




SNR analysis of TTR images

SNR = 6.18 (the average pixel intensity in the foreground over the average pixel intensity in the background)

Foreground pixels with original intensities






Background pixels with original intensities




frgr = pixels with intensity between 15 & 30
bkgr = pixels with intensity between 5 & 15

SNR analysis of CTL images







SNR analysis of CTL images


SNR = 2.24 (the average pixel intensity in the foreground over the average pixel intensity in the background)






frgr = pixels with intensity between 50 & 100
bkgr = pixels with intensity between 10 & 50

SNR analysis of CTL images





frgr = pixels with intensity between 50 & 100
bkgr = pixels with intensity between 10 & 50
SNR analysis of CTL images



frgr = pixels with intensity between 50 & 100

bkggr = pixels with intensity between 10 & 50

SNR analysis of CTL images

